

RB SYBR qPCR

Description:

The **RB SYBR qPCR Kit** is supplied as a 2X concentrated ready-to-use mixture for amplification of specific DNA fragments. It can be used in place of a single reagent for real-time quantitative PCR. The mixture contains Hot Start Taq DNA polymerase and other components required for real-time quantitative PCR with the exception of template and primers. Hot Start Taq DNA polymerase is complex mixture of a thermostable Taq DNA polymerase and specific monoclonal antibodies. It can eliminate amplification artifacts such as primer-dimer formation and mis-priming during preamplification stage and thus may provide improved specificity when compared to standard DNA polymerase. This not only saves valuable time in the laboratory, but also reduces the number of pipetting and reagent handling errors. A 1X RB SYBR qPCR mix contains 2 U of recombinant Hot Start Taq DNA polymerase, reaction buffer, 3 mM MgCl₂, 250 μM dNTP, SYBR Green I, ROX passive reference dye, and enzyme stabilizer sufficient to allow efficient amplification of template in a 25 μl reaction.

Components of the kit :

1 tube of 1 ml RB SYBR qPCR Mix.

General Procedure:

1. Prepare the reaction mix by combining the indicated volumes of components of **2XRB SYBR qPCR Mix**:

2.

Experimental Sample

| | |
|--------------------------------|---------|
| Nuclease-free H ₂ O | X μl |
| cDNA template | Y μl |
| Upstream primer (10 pM) | 1 μl |
| Downstream primer (10 pM) | 1 μl |
| 2XRB SYBR qPCR Mix | 12.5 μl |
| Total | 25 μl |

3. Place the tubes in a qPCR machine and proceed directly to amplify DNA. An initial denaturation step should be performed first at 94~95 °C for 5 minutes. This inactivates the antibody and initiates the hot-start PCR, then proceed a general PCR program.

Office, floor, noorbaran Building, 22Bahman St,Bozorgmehr Ave, Isfahan, Iran

www.rnabiotech.com

info@rnabiotech.com

sales@rnabiotech.com

03132673138-39

09336700380-390

Note:

1. Typical amplified sizes for real-time quantitative PCR products are less than 300 bp.
2. For either the two-step or the three-step protocol, set the cycle numbers in the range of 30~45 cycles.

**Examples of typical cycling parameters
Two-Step program**

| Step | Temp | Time | Detection |
|-------------------------|---------|----------|-----------|
| Denaturation | 94~95°C | 20~30s | off |
| Annealing/ Extension | TM-5°C | 30s~1min | on |

Three-Step program

| Step | Temp | Time | Detection |
|--------------|---------|------------|-----------|
| Denaturation | 94~95°C | 20 s~30 s | off |
| Annealing | TM-5°C | 30 s~45 s | off |
| Extension | 72°C | 30 s~1 min | on |